Uniblitz® FS6

6mm Uni-Stable Laser Shutter

Overview

The Uniblitz FS series shutters are designed and optimized to operate directly from +5VDC and do not require a separate driver. Removing the +5VDC (0VDC) closes these shutters. A simple control circuit can be used to operate these shutter devices from a TTL trigger pulse. This control can be also accomplished with our new **VLM1 TTL control interface module**, which is now available.

This low-cost innovation provides the reliability of Uniblitz shutters (typical lifetime >1 million cycles) at a single operating voltage.

Key Features

- 6mm aperture
- Default closed operation
- RoHS Compliant
- Transfer time on opening:3.3 milliseconds
- Transfer time on closing:6.0 milliseconds

Specifications


Electrical Specifications	
Coil resistance	12 OHMS
Voltage to Open	+5 VDC at 425 mA
Hold Voltage	+5 VDC

Weight Half-Housed 5 Weight Housed 9 Operating Temp.	5.7 g 9.5 g 4.5 g
Weight Housed 9 Operating Temp.	
Operating Temp.	4.5 g
operating remp.	8
Max Opening Bounce 1	40 - +65 C
max. opening bounce	5%
Max. Closing Bounce 5	%
Max. Freq. of Operation ¹ 5	00 Hz
Number of Shutter Blades 1	

¹ (Continuous/Burst) Continuous frequency rating specified at shutter's minimum exposure time (MET).



Shutter Timing Data



¹ Under no circumstances should any type of lubricant be applied to the shutter blade area. Lubricating the shutter blades will likely slow the shutter down and may eventually render it inoperable.

FS6 (w	1/ 5VDC and "T" blades) ¹ Time	(msec.)
O - A	Delay time on opening after current applied	5.2
A - C	Transfer time on opening	3.3
O - C	Total opening time	8.5
C - E	Min. dwell time with min. input pulse	3.7
B - F	Min. eqivalent exp. time	8.4
D - E	Delay time on closing after current removed	2.2
E - G	Transfer time on closing	6.0
A - G	Total window time	13.0
MET	Min. exposure time	10.0
TEP	Typical exposure pulse	12.0

Product Options

FS6 **23456**-7


1 Shutter Series:

Ex: FS6S2Z0L-EC

FS6: Normally Closed
FSR6: Normally Open
5 Electronic Sync:
C Licettonic Synta
0: Omitted
1: Included

2 Voltage:
S: Standard 5 VDC

3 Housing: 2 1: Un-Housed 2: Half-Housed 3: Fully-Housed



6 Connector: L: 18" Flying Leads

7 Encapsulated Coil: EC: Included 5 Leave blank if not required

² Various mounting methods available depending on housing option.


³ Other blade coating options may be available by special order.

⁴ Input side only; Teflon® coating is on opposite side to protect shutter blade surface. Light source must be input to the reflective side only.

⁵ Vacuum compatibility up to 10E-10 Torr.

Electronic Sync.

The synchronization system for FS shutter devices incorporates a small magnet mounted to the driving mechanism and a Hall effect sensor. When the device achieves approximately 80% of full open, the magnet causes the Hall effect sensor to change state, producing a signal to indicate that the shutter has switched to the active state. Shown to the left is the FS series shutter schematic which incorporates the electronic synchronization system. There is no connection to the designated synchronization pins when an electronic sync. is not selected.



Typical Control Circuit


This simple control circuit to the right can be used to operate the shutter device from a TTL trigger pulse. This control can also be accomplished with the <u>VLM1 interface</u> module.

¹ Note that the shutter will not operate directly from a TTL signal.

+5V. 0.5A 1. ALL RESISTORS 1/8W, 5%, EXCEPT AS NOTED. Red Wire 2. DIODES D1 & D2 PREVENT INPUT VOLTAGES EXCEEDING COIL (A) +5V AND DAMAGING MOSFET Q1. 3. +12V MAXIMUM INPUT FOR 1/8W RESISTORS. SHUTTER COIL DCR = 12 +5V Brown Wire COIL (B) D1 1N4148W 太 R1 Q1 TTL INPUT DMN2046U 1K OR EQUIVALENT D2 R2 本 1N4148W 10K


Uniblitz® FS6 Technical Drawings

Half-Housing - FS6S2T0L



- 1. UNITS: [MM] INCH.
- 2. LIMITED HIDDEN LINES SHOWN FOR CLARITY.
- 3. LEADS WIRES REMOVED FOR CLARITY.

Un-Housed - FS6S1T0L



Uniblitz® FS6 Technical Drawings

Full-Housing - FS6S3T0L

